

PROCARE SKIN BARRIER

Skin Barriers Composed of Skin-Friendly Natural Materials

Six types can be used for various purposes.



ÄLCARE

ALCARE Co., Ltd.
INTERNATIONAL SALES DIVISION

19F, Arcacentral Bldg.
1-2-1 Kinshi, Sumida-ku
Tokyo 130-0013, Japan
TEL:81-3-5611-7819
FAX:81-3-5611-7827
<http://www.alcare.co.jp/e>

Printed in China
Mar 2005

ÄLCARE


PROCARE SKIN BARRIER Series

Skin-friendly natural materials such as “Karaya Gum” and “Citrus Pectin” are used in the PROCARE SKIN BARRIER Series.

PROCARE SKIN BARRIER Series

PROCARE WAFER

Wafer Skin Barrier



- It is most suitable for partial use on wrinkles or hollows in the peristomal area.
- It can be freely cut to the shape of the hollow or the wrinkled area, or it can be used in layers if necessary.



Type	Product Code No.	Specification (Length×Width)	Contents Per Box
10	12581	100mm×100mm	10 pcs
20	12582	200mm×195mm	5 pcs

PROCARE POWDER

Powder Skin Barrier



- This powder fills the gap between the stoma and the face plate, protecting the skin from discharge.
- It is easy to apply, even to a small gap. The powder offers excellent water absorption.



Type	Product Code No.	Specification (Net Weight)	Contents Per Box
—	12611	50g	1 bottle

PROCARE PASTE

Paste Skin Barrier



- By applying this paste on the inside of the face plate, the adhesion of appliances can be improved.
- A laminate tube keeps the paste moist, making it easy to squeeze the tube down to the last.



Type	Product Code No.	Specification (Net Weight)	Contents Per Box
—	14431	50g	1 tube

PROCARE SOFT WAFER RING


Wafer Skin Barrier



- It is most suitable for partial use on wrinkles or hollows in the peristomal area.
- A moderately soft wafer that can be molded with the hands (no scissors needed).
- A ring-shaped wafer that is easy to apply to the peristomal skin area.
- A double-sided adhesive wafer that adheres well to the face plate of various appliances.



Type	Product Code No.	Specification (Diameter of Hole)	Contents Per Box	
Ring	20	15591	20mmφ	8 pcs
	30	15592	30mmφ	8 pcs



PROCARE MF PASTE

Paste Skin Barrier



- It fills the gap between the stoma and face plate, protecting the skin from discharge.
- An alcohol-free paste which will not make the skin sting but will remain soft.



Type	Product Code No.	Specification (Net Weight)	Contents Per Box
—	12601	100g	1 tube

PROCARE SOFT WAFER STICK

Wafer Skin Barrier



- It is most suitable for partial use on wrinkles or hollows in the peristomal area.
- An extremely soft wafer that can be freely molded.
- A double-sided adhesive wafer that adheres well to the face plate of various appliances.
- Ideal for application to uneven hard-to-apply surfaces such as the area around the fistula.



Type	Product Code No.	Specification (Length×Width)	Contents Per Box
Stick	15601	15mm×100mm	10 pcs

